Accuryn[®] Monitoring System





IAP ACCURACY | CAUTI RISK REDUCTION | WORKFLOW IMPROVEMENTS



The Accuryn[®] Solution provides vital information from the kidney that is accurate, continuous, and remote.

- Intra-abdominal pressure at the push of a button, with continuous trends
- Patented Accuryn SmartFoley[™] is designed to reduce the risk factors associated with CAUTI according to the CDC guidelines¹
- Continuous AKI surveillance²
- EMR Connectivity through Capsule[™]

Active Drain Technology: Air makes a u-turn at the foley barb and mechanically clears the drainage line. This patented feature addresses the need to "milk the foley" without applying negative pressure on the bladder.



IAP/UO/Temp

UO/Temp

Foley Adapter (6-24Fr)

EVEL

FLUID MANAGEMENT | AKI SURVEILLANCE | DATA-DRIVEN DECISIONS

The AccuTab[™] combined with the Accuryn Monitor, enables Providers to achieve their treatment goals of fluid balance while monitoring for early indications of renal failure.



Individualized care

The Accuryn system continually monitors urine output per the KDIGO guidelines (0.5ml/kg/hr). Simply enter the patient's weight and we take care of the rest.²

Acute kidney injury surveillance

AccuTab[™] monitors UO (in 15 minute increments) and Alerts if AKI criteria is reached per the KDIGO guidelines.



Annotate interventions

Clinicians can annotate interventions, and see near real time, the response of the kidney.

Furosemide Stress Test (FST)

AccuTab[™] monitors urine output against customized thresholds, and displays the results to enable timely interventions.

Operating room mode

Features a stop watch function to enable OR staff to start a UO counter during surgery.







Accuryn Monitor & AccuTab™

| Part Number | Description |
|---------------|--|
| FGN-06-102051 | Accuryn [®] Monitoring System |
| FGN-06-102052 | AccuTab [™] Monitoring System |
| FGN-06-102053 | Accuryn [®] Monitor and AccuTab™ System |

Accuryn SmartFoley[™] Catheters

| Part Number | Description |
|------------------|---|
| FGS-06-IAP.16FST | 16F Accuryn SmartFoley™ IAP/UO/T– 10 PK |
| FGS-06-UOT.16FST | 16F Accuryn SmartFoley™ UO/T – 10 PK |
| FGS-06-UOT-SMCS | Accuryn SmartFoley™ Adapter – 10 PK |
| FGN-06-100779 | Accuryn® Temperature Cable |

Accuryn CompleteCare Service

| Part Number | Description |
|---------------|------------------------------------|
| FGN-06-101693 | AccuTab™ Complete Care |
| FGN-06-102067 | Accuryn [®] Software Care |
| FGN-06-102068 | AccuTab [™] Software Care |

Accessories

| Part Number | Description |
|-----------------|--|
| FGN-06-100195 | Accuryn [®] Monitor Floor Stand |
| FGN-06-102111 | High Mount with Stand |
| FGN-06-102112 | High Mount |
| FGN-06-100202 | Accuryn [®] Monitor Mount with Clamp |
| FGN-06-100203 | Accuryn [®] Monitor Stand and Mount with Clamp |
| FGN-06-100800 | Accuryn® Monitor Mount w/ out Clamp |
| FGN-06-100823 | Accuryn [®] Monitor Stand and Mount w/out Clamp |
| FGN-06-3003 | Accuryn [®] Monitor 18V Power Supply |
| FGN-FLP-0008-17 | AccuTab [™] Swivel Mount |
| FGN-PRC-0003-08 | AccuTab [™] Flat Mount |
| FGN-06-101679 | AccuTab™ USB Cable |
| FGN-06-101678 | AccuTab™ Power Cable |
| FGN-MD100DS | AccuTab [™] Docking Station |
| FGS-06-101623 | Centurion Securement System |
| FGS-06-101796 | Grip-Lok Catheter Securement Device, Sterile |

References:

- 1. CDC Guidelines to reduce CAUTI
- https://www.cdc.gov/infectioncontrol/guidelines/cauti/
- 2. Kidney International Supplements (2012) 2, 4; doi:10.1038/kisup.2012.4

24/7 Product helpline 1.833. Accuryn 1.833.222.8796

www.potreromed.com info@potreromed.com 1-888-635-7280

Potrero

MKT-06-102066 Rev D